

## Let's Enjoy

Shower Chair

Article no. LEN101 Black

Article no. LEN201 White

Let's Enjoy is a new shower chair in a modern design intended for use in the shower or the bathroom. The chair features a beautiful and stable construction for safe and secure use. The legs are easy to adjust to desired height with one easy motion. The chair is intended for intimate cleaning.



Let's Enjoy is adjustable from 42 to 53 cm. The height is adjusted by pressing the buttons on the legs and pulling the lower legs to the desired positions, after which the buttons will be locked in place in the adjustment holes.

**Important:** Always adjust the legs to the same height to avoid falling accidents.

**Technical data:**

Material:

Seat ABS

Tube fram aluminum

Feat TPE

Height adjustable from 42 to 53 cm


Weight 3 kg

Tested according to ISO EN 12182:2012


Max user weight 150 kg

**Warranty:**

Two years warranty



42 - 53 cm



52 cm

37 cm