

Let's Shower

Shower Stool

Article no. LSH101 Black

Article no. LSH201 White

Let's Shower is a new shower stool in a modern design intended for use in the shower or the bathroom. The stool features a beautiful and stable construction for safe and secure use. The legs are easy to adjust to desired height with one easy motion. The stool is easy to fit in any space thanks to its small shape and round seat.



Let's Shower is adjustable from 45 to 53 cm. The height is adjusted by pressing the buttons on the legs and pulling the lower legs to the desired positions, after which the buttons will be locked in place in the adjustment holes.

Important: Always adjust the legs to the same height to avoid falling accidents.

Technical data:

Material:
Seat ABS
Tube fram aluminum
Feat TPE
Height adjustable from 45 to 53 cm
Weight 2 kg
Tested according to ISO EN 12182:2012
Max user weight 150 kg

Warranty:
Two years warranty

