

Let's Step

Foot Stool

Article no. LST101

Let's Step is a new foot stool with a modern design intended for use in the shower or the bathroom. The stool has a beautiful and stable construction for safe and secure use.

The stool has ribs on the top to give a secure grip when the stool is wet and it is black to give high visibility against bathtubs.



Technical data:

Material:

Seat ABS

Tube frame aluminum

Feet TPE

Weight 2,8 kg

Tested according to

ISO EN 12182:2012

Max user weight 150 kg

Warranty:

Two years warranty

