

## Let's Go Out

Outdoor rollator

Black Article No. 0120

Red Article No. 0130

Beige Article No. 0150



Let's Go Out rollator is light and steady. It weighs only 6.1 kgs and has a width of only 58.5 cms. It is easy to fold and the bag is folded together with the rollator. Let's Go Out is available in three colours: Black/silver with seat in black artificial leather, Beige/silver or Red with seat in brown artificial leather. The bags also match the rollator in colour. The large front wheels make it easy to turn and climb kerbs. Let's Go Out has fenders on the rear wheels and built-in spring suspension that provides comfortable walking. The rubber handles are ergonomically shaped and can be adjusted in height from 78 to 96 cms to fit different users.

Let's Go Out has a spring suspension in the rear wheels for smoother walk.



When Let's Go Out is folded it is easy to bring with you in the car or on a bus as the width is only 20 cm.



Bag. Brown detail Art. No 151  
Black detail Art. No 121



Black plastic tray. Art. No 0381



Black or Brown backrest band. Art. No 0383



Crutch holder. Easy to mount without any tools, only snap on. Art. No 0382



### Technical specifications:

Tested according to ISO 11199-2:2005 for a maximum user weight of 130 kg.  
Width 58.5 cm.  
Length 67 cm.  
Weight 6.1 kg.  
Seat height 59.5 cm.

Front wheel diameter 24.5 cm.  
Rear wheel diameter 19.5 cm.  
Rear wheels with spring suspension.  
Handle height adjustable from 78 to 96 cm.  
Maximum bag load 10 kg.

