



Let's Shop. Spare Parts.

Article No. LS101 Black
Article No. LS102 Grey

06.2019




Rubber handle with tube, reflector and end stop.




Art. No. LS307 Left.
Art. No. LS308 Right.

Brake levers, adjustment screws and brake housing.




Art. No. LS306.
Set Right + Left.

Bushings with bottom tube stoppers, set.




Art. No. LS320.
Set Right + Left.

Height adjustment buttons, set.




Art. No. LS319.
Set Right + Left.

Mounting set for seat. 2 pcs.




Art. No. LS324.
Set Right + Left.

Dome nuts for seat mounting. Set 4 pcs.




Art. No. LS329.
Set 4 pcs.

Seat.




Art. No. LS310.

Unfold hook with mounting screws and spring.




Art. No. LS311.

Linkage arm with reinforcement, screws and hook.




Art. No. LS312 Left.

Linkage arm with reinforcement and screws.




Art. No. LS321 Right.

Middle screw for folding cross, complete.




Art. No. LS326.

Brake wire guide inside rear leg. Set 2 pcs.




Art. No. LS313.
Set Right + Left.

Front fork with front fork bearings, set 2 pcs.




Art. No. LS318.
Set Right + Left.

Front fork bearings, set 2 pcs.




Art. No. LS322.
Set Right + Left.

Rear wheel bracket with screws.




Art. No. LS323.
Set Right + Left.

Brake yokes.




Art. No. LS328.
Set Right + Left.

Wheel with bearings and mounting set.




Art. No. LS317.

Mounting set for front and rear wheels. Set 4 pcs.




Art. No. LS314.
Set Front + Rear.

Ball bearings for front and rear wheels. Set 4 pcs.



Art. No. LS477.
Set 4 pcs.

Brake wire with spring for left and right side.



Art. No. LS305.
Set Right + Left.