

# Let's Help Me


Commode

Rose Gold Art. No. LHM100

Silver Art. No. LHM110



A commode helps you get up smoothly from the toilet seat. With Let's help Me, you can set the seat height to your needs. The telescopic legs and anti-slip rubber feet ensure stability on any surface. Its ergonomic back and armrests aid in getting up and down. As this product is also water-resistant, it can also be used as a stool in the bathroom or shower.



### Technical specifications:

Tested according to: ISO 17966:2016, EN ISO 10993-1 for a maximum user weight of 160 kg.

Height	64 - 81 cm
Seat height	41 - 58 cm
Depth	51 - 59 cm
Width	59 - 64 cm
Weight	3.1 kg



### Box measurements

57,8 cm x 19,7 cm x 57,8 cm  
4,53 kg

### Container size

20 GP	420 sets
40 GP	892 sets
20 HQ	952 sets

