

Let's Sing

Shower chair

White Art. No. LSG101

Ice Blue Art. No. LSG301

Copper Art. No. LSG401



Let's Sing is a new shower stool in a modern design intended for use in the shower or the bathroom. The stool features a beautiful and stable construction for safe and secure use. The legs are easy to adjust to desired height with one easy motion.

A wire basket is available as an accessory that should be placed under the seat for storing shower articles or other products.

Let's Sing is adjustable from 45 to 56 cm in height. The height can be adjusted by pressing the buttons on the legs and then moving the lower tubes to the preferred height, the buttons will then pop back into place and lock the legs.

Important: Always adjust all the legs to the same height to avoid falling accidents.



Ice Blue. Art. No. LSG301



Copper Art. No. LSG401

Technical specifications:

Tested according to ISO 12182:2012

Max. user weight 150 kg

Adjustable height 45 – 56 cm

Weight 2.4 kg

Material:

Seat ABS

Tube frame Aluminum

Feet TPE

