

Let's Get Up

Safety bed handle

White Art. No. LGU101



Let's Get Up provides protection against accidental falls from the bed and offers you dependable assistance in transitioning to and from the bed.

The ergonomic handle is covered with soft rubber and can be adjusted in various heights. Let's Get Up is equipped with an adjustable bed frame bracket and fits most traditional home and medical beds.

This product is a fantastic option for people who lack balance or simply need a little extra assistance when getting into bed.

Maximum user weight: 150 kg.

Technical specifications:

Tested according to IEC 60601-2-52:2009 for a maximum user weight of 150 kg.

| | |
|--------------------|-------------|
| Handle height | 81 – 97 cm. |
| Bed bracket height | 9 – 45 cm. |
| Floor plate length | 59 cm. |
| Floor plate width | 35 cm. |
| Weight | 11.9 kg. |

