

Let's Enjoy

Shower chair

Black Art. No. LEN101
 White Art. No. LEN201
 Ice Blue Art. No. LEN301
 Copper Art. No. LEN401



Let's Enjoy is a new shower chair in a modern design intended for use in the shower or the bathroom. The chair features a beautiful and stable construction for safe and secure use. The legs are easy to adjust to desired height with one easy motion. The chair is intended for intimate cleaning.

Let's Enjoy is adjustable from 42 to 53 cm. The height is adjusted by pressing the buttons on the legs and pulling the lower legs to the desired positions, after which the buttons will be locked in place in the adjustment holes.

Important: Always adjust all the legs to the same height to avoid falling accidents.



White. Art. No. LEN201



Ice Blue. Art. No. LEN301



Copper. Art. No. LEN401

Technical specifications:

Tested according to ISO 12182:2012
 Max. user weight 150 kg
 Adjustable height 42 – 53 cm
 Weight 2.4 kg



Material:
 Seat ABS
 Tube frame Aluminum
 Feet TPE

