

Let's Move

Outdoor rollator

Black Art. No. LM201

Grey Art. No. LM202



Let's Move, stylish and innovative outdoor rollator from Trust Care. The rollator has a sturdy lightweight aluminum frame and whole rollator weighs only 5.4 kg. With such low weight it is easy and pleasurable to use the rollator in any conditions. Of course it has an effective brake and parking brake that are easy to operate. A fabric bag is also available for Let's Move. The rollator is easy to fold and can stand by itself when folded. Folded width is only 20 cm, such feature makes it easy to store or take it with you wherever you go. Let's Move has ergonomically designed rubber handles, adjustable from 78 to 93 cm in order to fit different users.

Let's Move is tested and approved according to ISO specifications for 130 kg user weight.



Black. Art. No. LM201



Grey. Art. No. LM202

Technical specifications:

Tested according to ISO 11199-2:2005

Max. user weight	130 kg
Width	58 cm
Length	65 cm
Weight	5.4 kg
Handle height is adjustable from	78 – 93 cm
Seat height	60 cm
Wheel diameter	20 cm
Max. bag load	12 kg

