

# Let's Sit Up

Armchair

|       |                 |                          |
|-------|-----------------|--------------------------|
| Grey  | Art. No. LSU201 | EAN code: 07350036885090 |
| Green | Art. No. LSU301 | EAN code: 07350036885106 |
| Beige | Art. No. LSU401 | EAN code: 07350036885113 |

Let's Sit Up is a chair inspired by classic Scandinavian design. Its main function is the ability to help the user stand up with the help of a built in mechanism under the seat. The seat raises upwards with the push of a lever and slowly raises the user up. The mechanism is powered by a gas spring and there is no need to motors or electricity. Its easy to place anywhere in the room without having to worry about cables across the floor. The compact design of the mechanism and the shape of the chair makes it look like a regular armchair.

The fabric is of high durability, flame-proof and you can choose from different colors to fit every household. The chair is made of moulded foam and aluminum that has been coated with wood grain film to make it look and feel just like wood. The built in handles gives extra support to the user when standing up.



**Technical data:**

|                   |        |
|-------------------|--------|
| Max. user weight: | 110 kg |
| Width:            | 58 cm  |
| Length:           | 74 cm  |
| Height:           | 89 cm  |
| Weight:           | 25 kg  |
| Seat height:      | 46 cm  |
| Seat width:       | 50 cm  |
| Seat depth:       | 48 cm  |

