

# Let's Frisbee

Shower chair

Black Art. No. LFR101  
Ice Blue Art. No. LFR301  
Copper Art. No. LFR401



Let's Frisbee is a new bath stool with modern design, easy to use in the shower or in the bath tub. The stool has a beautiful and unique construction for safe and stable use. It also has a rotatable soft seat which makes it easy for the user to handle their hygiene in a comfortable manner. The stool is designed to fit in a variety of places with its small and round shape.

Let's Frisbee is adjustable from 45 to 53 cm in height. The height can be adjusted by pressing the buttons on the legs and then moving the lower tubes to the preferred height, the buttons will then pop back into place and lock the legs.

**Important:** Always adjust all the legs to the same height to avoid falling accidents.



Ice Blue. Art. No. LFR301



Copper. Art. No. LFR401

**Technical specifications:**

Tested according to ISO 12182:2012  
Max. user weight 150 kg  
Adjustable height 45 – 53 cm  
Weight 2.5 kg



Material:  
Seat ABS  
Tube frame Aluminum  
Feet TPE

