

Let's Twist Again

Exploded View and Part Numbers

- LTA300B Top tube cap Black
- LTA300S Top tube cap Silver
- LTA301Y Elbow arm support Yellow
- LTA301O Elbow arm support Orange
- LTA301G Elbow arm support Green
- LTA301T Elbow arm support Transparent
- LTA301B Elbow arm support Black
- LTA302B Top screw Black
- LTA302S Top screw Silver
- LTA303B Crutch body tube Black
- LTA303S Crutch body tube Silver
- LTA304 Handle inner tube
- LTA305Y Handle left Yellow
- LTA305O Handle left Orange
- LTA305G Handle left Green
- LTA305T Handle left Transparent
- LTA305B Handle left Black
- LTA306Y Handle right Yellow
- LTA306O Handle right Orange
- LTA306G Handle right Green
- LTA306T Handle right Transparent
- LTA306B Handle right Black
- LTA307 Handle reinforcement rod
- LTA308 Handle screw
- LTA309Y Handle cap Yellow
- LTA309O Handle cap Orange
- LTA309G Handle cap Green
- LTA309T Handle cap Transparent
- LTA309B Handle cap Black

- LTA310 Adjustment button holder
- LTA311 Adjustment tube cap
- LTA312 Adjustment lock collar
- LTA313B Adjustment lock ring Black
- LTA313S Adjustment lock ring Silver
- LTA314B Adjustment tube Black
- LTA314S Adjustment tube Silver
- LTA315 Foot

